CODE: FIT 132 FIT 176

Food Inflammation Test (FIT) KBMO DIAGNOSTICS

THE PROCESS

- The blood sample can be collected on any day
- This is not a fasting test. Eat your usual diet before the test
- Read the instructions in the collection kit
- KBMO 132 has one test card (5 spots) and KBMO 176 has two cards (10 spots)

BLOOD FLOW TIPS

- 10-20 minutes prior to taking your sample:
 - o drink 2-4 glasses of water or tea
 - o ensure your hands are warm
 - o move your body to get your blood pumping
- Warm your hands by rubbing them together and shaking them
- Rub the finger you will use to stimulate blood flow
- Having your hand below your heart can help with blood flow
 See next page for more tips

4 Weeks before the test stop taking:

• Immunosupressive drugs (including topical steroids)

Day before

• ensure you are well hydrated



READ THE INSTRUCTIONS INCLUDED IN THE TEST KIT BEFORE COMPLETING TEST

Ensure sample is collected correctly.

The lab charges a \$50 retesting fee for incorrect samples





CODE: FIT 132 FIT 176

Food Inflammation Test (FIT) KBMO DIAGNOSTICS

COLLECTION TIPS

- Use the alcohol swap to disinfect your finger before and straight after the lancet prick, then start collection
- Do not touch the test circles this will contaminate the sample and result in a re-test fee
- Do not place more blood on the card once dry
- You must fully fill each circle on the collection card







Not all circles to the edge of dotted lines

POSTING

- You will need to send your samples to the lab using the postage satchel provided
- Samples must be completely air dried before posting
- Post within 1-2 days of completing
- Notify your health coach via slack when you have posted share a photo of the shipping label on the postage packaging, this will allow us to track the sample if required



READ THE INSTRUCTIONS INCLUDED IN THE TEST KIT BEFORE COMPLETING TEST

Ensure sample is collected correctly.

The lab charges a \$50 retesting fee for incorrect samples

