

CODE: DUN009

Precision sIgA

SALIVA ANTIGEN TEST

THE PROCESS

- The saliva sample must be **collected on Monday**
- You will need to send your samples to RN Labs using the postage satchel provided.
- Send samples on **Monday or Tuesday only**
- Samples must remain **refrigerated until ready to send**.
- **Do not freeze samples**. Frozen samples will be rejected.
- **Notify your health coach** when you have completed the test via Slack - **share a photo of the shipping label** on the postage packaging, this will allow us to track the sample if required

4 Weeks before the test stop taking:

- **Immunosuppressive drugs (including topical steroids)**

10 minutes before test,

- do not eat or drink



READ THE INSTRUCTIONS INCLUDED IN THE TEST KIT BEFORE COMPLETING TEST

Ensure sample is collected correctly.
The lab charges a \$50 retesting fee for incorrect samples